


FRONT VIEW



BACK VIEW

1. P/N 130-BIS-0007.
2. GENERAL SPECIFICATIONS:
 - 2.1 VISUAL INSPECTION:
 - 2.1.1. NO RUPTURES OR DAMAGE.
 - 2.1.2. NO MATERIAL CUTS.
 - 2.1.3. NO CUTTING NICKS.
 - 2.1.4. NO INCORRECT MARGINS.
 - 2.1.5. NO COLORATIONS OR STAINS AS PER PANTONE 2612 C.
 - 2.1.6. NO ILLEGIBLE WORDS.
 - 2.1.7. WRITING SURFACE MATERIAL PANEL.


2.1.8. WRITING SURFACE MATERIAL PANEL.

■ PANTONE 2612 C. (Glossy texture)

3 DIMENSION INSPECTION

4. FOR ESTABLISHMENT LABS ONLY:

4.1 INSPECT PER QI-011.

Designed by M.VALVERDE	Release Date Refer to CQ	Tolerance/ unless specified TWO PLACE DECIMAL ±1,00 mm	MOTIVA Patient Identification Card with Qid.		
 Establishment Labs Property - All Rights Reserved		Released per: CHG-000464	DRW No. DRW-544	Revision 4	Sheet 1 / 1
		Dimensions mm unless specified	Scale 1:1		