

Magnetic Resonance Imaging (MRI) Information

LVIS[®] EVO[™]

Intraluminal Support Device

*Please carry your card
at all times and show it
to any medical personnel
who may be treating you.*

The LVIS[®] EVO[™] device has been shown to be MR conditional in MRI systems operating under the following conditions: static magnetic field of 3.0 Tesla or less.

Before you have an MRI scan, or for questions regarding your LVIS[®] EVO[™] device implant, please contact your physician or MicroVention at 1-800-990-8368.

Patient Name _____

Telephone _____

Date of Implant _____

No. of Stents Implanted _____

Site of Implant(s) _____

Implanting Physician _____

Hospital _____

Address _____

Telephone _____

Primary Physician _____

Telephone _____

MICROVENTION is a registered trademark of MicroVention, Inc. LVIS® EVO™ is a trademark of MicroVention, Inc. Scientific and clinical data related to this document are on file at MicroVention, Inc. Refer to Instructions for Use for additional information. Federal (USA) law restricts this device for sale by or on the order of a physician. ©2019 MicroVention, Inc. IFU100053 Rev A 07/19

LVIS® EVO Device (affix label if available):
Model Number:
Size:
Lot Number:

Patient Implant Card



LVIS® EVO™

Intraluminal Support Device

1311 Valencia Avenue, Tustin, CA 92780 USA

1-800-990-8368 • www.microvention.com